

Liposuction Cannulas

EFFICIENT BODY SHAPING WITH PRECISE DISPOSABLE CANNULAS

Classics Cannula Type with anatomical handle
in high-quality single-use version

Major 5 large openings, one closed side

Speed+ 3 large and 12 smaller openings

Pyramid 3 large openings, one closed side

Inline 2 large openings, only on one side

Harvester 24 smaller openings

Different Lengths and Thicknesses can be available

Tumescent Infiltration Cannulas
with Luer Lock connection

Also available
with threaded connection
for LipoSurg system and
manual Handles

Liposuction Cannulas with Anatomical Handle

Typ	REF	Ø	x Length
Major	LX5P45300	4.5 mm	x 300 mm
Speed+	LXS+40300	4.0 mm	x 300 mm
Pyramid	LX3P40300	4.0 mm	x 300 mm
	LX3P30250	3.0 mm	x 250 mm
Inline	LX2L30250	3.0 mm	x 250 mm
	LX2L30130	3.0 mm	x 130 mm
Harvester	LXH40300	4.0 mm	x 300 mm
	LXH30250	3.0 mm	x 250 mm

Tumescent Infiltration Cannula, Luer Lock

LXT20180	2.0 mm x 180 mm
LXT20250	2.0 mm x 250 mm
LXT25300	2.5 mm x 300 mm

Lipo Transfer Cannulas



WITH LUER LOCK CONNECTOR

LIPO HARVESTER


for evenly structured Lipoaspirate in manual extraction

SINLE PORT CANNULA

for the precise, atraumatic Deposition in the target tissue

V-DISSECTOR

for the resolution of scarred indents and fibrous tissue



Lipo Harvester

Typ	REF	x Length
Micro	LXHmi20150	2.0 mm x 150 mm
Curette	LXHS25200	2.5 mm x 200 mm
Atraumatic	LXHA25200	2.5 mm x 200 mm
Multiport	LXH30250	3.0 mm x 250 mm

SINGLE PORT CANNULA

Typ	REF	x Length
LX Spoon	LXiS20150	2.0 mm x 150 mm
LX	LXi20150	2.0 mm x 150 mm
	LXi25180	2.5 mm x 180 mm
V-Diss	LXV15090	1.5 mm x 90 mm
	LXV20150	2.0 mm x 150 mm

