




+ SOLUTIONS: MANUFACTURED IN SWITZERLAND +

Single-Use Cannulas

+ NOUVAG AG

St. Gallerstrasse 25
9403 Goldach
Switzerland
T +41 (0) 71 846 66 00
info@nouvag.com

+ NOUVAG GmbH

Schulthaißstrasse 15
78462 Konstanz
Germany
T +49 (0) 7531 1290-0
F +49 (0) 7531 1290-12
info-de@nouvag.com

Aspiration Cannulas



SINGLE PORT

REF		4643	4644
Cannula Ø	mm	2.0	3.0
Length	mm	150	300



THREE PORT

REF		4647	4648	4649
Cannula Ø	mm	2.0	3.0	4.0
Length	mm	200	300	300



MERCEDES

REF		4653	4654	4655	4656
Cannula Ø	mm	2.0	3.0	3.0	4.0
Length	mm	150	250	300	250



TWO PORT

REF		4645	4646
Cannula Ø	mm	2.5	3.0
Length	mm	200	300



ONE CENTRAL TWO LATERAL

REF		4650	4651	4652
Cannula Ø	mm	3.0	4.0	5.0
Length	mm	300	300	300



DOUBLE MERCEDES

REF		4657	4658
Cannula Ø	mm	4.0	5.0
Length	mm	250	300



TRIPPLE MERCEDES

REF	4659
Cannula Ø	mm 5.0
Length	mm 300



DOUBLE BASKET

REF	4662	4663	4664
Cannula Ø	mm 4.0	mm 5.0	mm 5.0
Length	mm 300	mm 300	mm 250



SATTLER

REF	4666	4667	4668	4669	4670
Cannula Ø	mm 2.0	mm 3.0	mm 3.0	mm 3.0	mm 4.0
Length	mm 150	mm 200	mm 300	mm 300	mm 300
Hole	30	30	30	30	30



SINGLE BASKET

REF	4660	4661
Cannula Ø	mm 3.0	mm 4.0
Length	mm 250	mm 300



CANDY CANE

REF	4665
Cannula Ø	mm 4.0
Length	mm 300



TOLEDO V

REF	4671	4672
Cannula Ø	mm 2.0	mm 3.0
Length	mm 150	mm 250

Infiltration Cannulas



BAYONETTE

REF		4673	4674
Cannula Ø	mm	3.0	4.0
Length	mm	300	300
Hole		10	10



ONE HOLE SPATULA

REF		4675
Cannula Ø	mm	2.0
Length	mm	150



30 HOLE

REF		4676	4677	4678	4679
Cannula Ø	mm	2.0	2.0	2.0	3.0
Length	mm	150	250	300	250
Hole		30	30	30	30